# Synthesis 1

**Betriebsanleitung**  
**Instruction Manual**  
**Mode d'emploi**  
**Instrucciones de empleo**  
**Istruzioni per l'uso**

<b>D</b>	DEUTSCH	Seite	3 - 55
<b>E</b>	ENGLISH	page	56 - 106
<b>F</b>	FRANCAISE	page	107 - 159
<b>ES</b>	ESPAÑOL	página	160 - 213
<b>I</b>	ITALIANO	pagina	214 - 266

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Niederspannungs-Richtlinie (73/23/EWG):

EN 61010

Thank you for buying a HEIDOLPH product. This appliance has been designed and tested in compliance with DIN EN ISO 61010. It features ease of handling and long service life.

## SUMMARY

<b>SUMMARY</b> .....	<b>56</b>
<b>STANDARD EQUIPMENT AND OPTIONS</b> .....	<b>58</b>
<b>GENERAL INFORMATION</b> .....	<b>60</b>
<b>SAFETY INFORMATION</b> .....	<b>61</b>
<b>SET UP</b> .....	<b>62</b>
1. Liquid-Phase Synthesis .....	62
1.1 Reaction Vessels .....	62
1.1.1 Setup.....	62
1.1.2 Install head.....	63
1.1.3 Built-in rotary valve.....	63
1.1.4 Septa.....	64
1.1.5 Capacity of reaction vessels.....	64
1.2 Locate reaction vessels in heating- and condensation block .....	64
1.2.1 Remove and install reaction vessels .....	64
1.2.2 Vacuum-/ protective gas fitting with distributor (plenum).....	67
1.3 Connect vacuum- / inert gas lines .....	69
1.3.1 Connect vacuum- / inert gas lines .....	69
1.3.2 3-way valve, vacuum- / inert gas (standard equipment).....	69
1.3.3 Valve settings, vacuum-/ inert gas .....	70
1.4 Less reaction vessels .....	70
2. Solid-phase synthesis .....	71
2.1 Reaction vessels .....	71
2.1.1 Setup:.....	71
2.1.2 Install and disassemble head / base .....	72
2.1.3 Base with built-in filter .....	72
2.1.4 Septa.....	72
2.1.5 Tube fitting, head and base.....	73
2.1.6 Open and close valve of reaction vessel .....	73
2.1.7 Capacity of reaction vessels.....	74
2.2 Locate reaction vessels in heating- and condensation block .....	74
2.2.1 Install vessels.....	74
2.3 Connect vacuum / inert gas with vessels / plenum.....	75
2.4 Connect vacuum- / inert gas to plenum ports.....	75
2.4.1 3-way valve, vacuum-/ inert gas (standard equipment).....	76
2.4.2 Valve settings, vacuum- / inert gas .....	76
2.5 Connect waste cube, solid-phase synthesis.....	77
2.5.1 Convert waste cube for deviating number of vessels .....	78
2.5.2 Schematic, arrangement of waste cube tubes .....	79
2.6 Less reaction vessels .....	79
3. Coolant fitting, liquid- or solid-phase synthesis .....	80
3.1 Coolant connector, condensation zone .....	80
3.2 Coolant fitting, conditioning zone for reactions below room temperature.....	80
4. Multiple-Evaporator.....	80

4.1	Reaction vessels .....	80
4.2	Locate reaction vessels in heater block .....	80
4.3	Connect vacuum .....	81
4.4	Set-up Synthesis 1 glassware kit for Multiple-Evaporator .....	81
4.4.1	Set-up stand.....	82
4.4.2	Install cooler .....	82
4.4.3	Connect vacuum and coolant.....	83
4.4.4	Install receiving flask .....	83
4.5	Schematic, Multiple-Evaporator without vacuum control .....	84
4.6	Schematic, Multi-Evaporator with vacuum control.....	84
5.	Electric wiring.....	86
6.	Install guard hood .....	86
6.1	Guard hood .....	86
6.2	Installation .....	87
7.	Exchange reaction blocks .....	87
	<b>OPERATION .....</b>	<b>89</b>
8.	General .....	89
9.	Controls, brief description .....	89
9.1	Power switch .....	89
9.2	Start / Stop shaking action .....	89
9.3	Start / Stop heater .....	90
9.4	Set / edit parameters .....	90
9.5	Set or display .....	90
9.6	Select value.....	90
9.7	Set temperature gradient .....	90
9.8	Activate temperature gradients .....	90
9.9	Identical .....	90
10.	Setting work in detail.....	91
10.1	Set shaking frequency.....	91
10.2	Set temperature (4 heater zones with different temperatures) .....	91
10.2.1	Setting.....	91
10.2.2	Temperature setting, limiting values.....	91
10.3	Set temperature (4 heater zones with same temperature) .....	91
10.4	Set timer.....	91
10.5	Set temperature gradient .....	92
10.5.1	Set values .....	92
10.5.2	Delete values .....	93
10.5.3	Limiting temperatures.....	93
10.6	Activate temperature gradient .....	93
10.7	Turn shaking action ON .....	94
10.8	Turn heater ON .....	94
11.	External temperature sensor.....	94
11.1	Temperatures below solvent's boil point .....	95
11.2	Temperatures above solvent's boil point.....	95
11.3	In the use as multi evaporator .....	95
12.	Setup (timer and temperature sensor) .....	96
13.	Filter action, solid-phase synthesis .....	97
14.	Use condensation zone .....	98
15.	Use heating zones with recirculating cooler (low temperature).....	98
16.	Access to rear reaction vessels .....	99
17.	Interface RS 232 .....	100

**CLEANING & SERVICING.....102**  
**UNINSTALL, FORWARD & STORE.....102**  
**DISPOSAL .....103**  
**TROUBLESHOOTING .....103**  
**SPECIFICATIONS .....104**  
**WARRANTY, LIABILITY & COPYRIGHT.....105**  
**QUESTIONS / REPAIR WORK.....105**

**important note**

**power supply (mains)**

**Caution**

**Caution: fire- or explosion hazard**

**Caution: danger of burns**

**repair & service information**

## STANDARD EQUIPMENT AND OPTIONS

	item	qty.	P/N 230V / 50/60Hz	P/N 115V / 50/60Hz
	Synthesis 1 Solid 16	1	552-14200-00	552-14200-01
or	Synthesis 1 Solid 20	1	552-12500-00	552-12500-01
or	Synthesis 1 Solid 24	1	552-10800-00	552-10800-01
or	Synthesis 1 Liquid 12	1	552-25000-00	552-25000-01
or	Synthesis 1 Liquid 16	1	552-22500-00	552-22500-01
or	Synthesis 1 Liquid 24	1	552-21000-00	552-21000-01
or	Synthesis 1 Multiple Evaporator 12	1	552-35000-00	552-35000-01
or	Synthesis 1 Multiple Evaporator 16	1	552-32500-00	552-32500-01
or	Synthesis 1 Multiple Evaporator 24	1	552-31000-00	552-31000-01
or	Synthesis 1 Solid System 16	1	552-64200-00	552-64200-01
or	Synthesis 1 Solid System 20	1	552-62500-00	552-62500-01
or	Synthesis 1 Solid Systems 24	1	552-60800-00	552-60800-01
or	Synthesis 1 Liquid System 12	1	552-85000-00	552-85000-01
or	Synthesis 1 Liquid System 16	1	552-82500-00	552-82500-01

or	Synthesis 1 Liquid System 24	1	552-81000-00	552-81000-01
or	Synthesis 1 Multiple Evaporator System 12	1	552-75000-00	552-75000-01
or	Synthesis 1 Multiple Evaporator System 16	1	552-72500-00	552-72500-01
or	Synthesis 1 Multiple Evaporator System 24	1	552-71000-00	552-71000-01
	Instruction Manual	1	01-005-004-41	01-005-004-41
	power cord	1	14-007-003-81	14-007-003-89

### Accessories

	item	P/N
	Reaction Block Solid 16	552-14210-00
	Reaction Block Solid 20	552-12510-00
	Reaction Block Solid 24	552-10810-00
	Reaction Vessel Solid 16	552-14230-00
	Reaction Vessel Solid 20	552-12530-00
	Reaction Vessel Solid 24	552-10830-00
	Filter 16	23-30-01-03-68
	Filter 20	23-30-01-03-69
	Filter 24	23-30-01-03-70
	Septa Solid	23-30-01-03-67
	PTFE-Tube	23-30-01-03-66
	Clamp Ring for PTFE-Tubes	23-30-01-03-71
	Temperature Sensor Probe (135mm) for Reaction Block Solid 24 / 20 / 16 and Liquid 24	552-50090-00
	Temperature Sensor Probe (170mm) for Reaction Block Liquid 16 and 12	552-50092-00
	Guard Hood	552-50091-00
	Waste Cube 16	552-44200-00
	Waste Cube 20	552-42500-00
	Waste Cube 24	552-40800-00
	Insert for Waste Cube 16	552-44210-00
	Insert for Waste Cube 20	552-42510-00
	Insert for Waste Cube 24	552-40810-00
	Test Tube Rack 16 for Waste Cube 16	552-44220-00
	Test Tube Rack 20 for Waste Cube 20	552-42520-00
	Test Tube Rack 24 for Waste Cube 24	552-40820-00
	Test Tube for Test Tube Rack 16	552-44260-00
	Test Tube for Test Tube Rack 20	552-42560-00

	<b>Test Tube for Test Tube Rack 24</b>	<b>552-40860-00</b>
	<b>Waste Cube Vessel</b>	<b>552-40060-00</b>
	<b>Blanking Screw Cap</b>	<b>23-30-01-03-72</b>
	<b>Reaction Block Liquid 12</b>	<b>552-25010-00</b>
	<b>Reaction Block Liquid 16</b>	<b>552-22510-00</b>
	<b>Reaction Block Liquid 24</b>	<b>552-21010-00</b>
	<b>Reaction Vessel Liquid 12</b>	<b>552-25030-00</b>
	<b>Reaction Vessel Liquid 16</b>	<b>552-22530-00</b>
	<b>Reaction Vessel Liquid 24</b>	<b>552-21030-00</b>
	<b>Glass Vessel Liquid 12</b>	<b>552-25060-00</b>
	<b>Glass Vessel Liquid 16</b>	<b>552-22560-00</b>
	<b>Glass Vessel Liquid 24</b>	<b>552-21060-00</b>
	<b>Septa Liquid</b>	<b>23-30-01-03-73</b>
	<b>Glassware set Synthesis 1</b>	<b>552-50070-00</b>
	<b>PTFE-Vacuum Tube</b>	<b>23-30-01-03-15</b>
	<b>PTFE-Connecting Tube 12 &amp; 16</b>	<b>552-28680-00</b>
	<b>PTFE- Connecting Tube 24</b>	<b>552-26480-00</b>
	<b>VAC control automatic RS 230V 50/60Hz</b>	<b>591-00341-00</b>
	<b>VAC control automatic RS 100/115V 50/60Hz</b>	<b>591-00341-01</b>
	<b>Vacuum valve</b>	<b>591-24000-00</b>
	<b>Woulff bottle</b>	<b>591-00071-00</b>
	<b>Tool Liquid 12</b>	<b>04-06-01-21-07</b>
	<b>Tool Liquid 16</b>	<b>04-06-01-21-08</b>
	<b>Tool Liquid 24</b>	<b>04-06-01-21-09</b>
	<b>Tool gas tube connector 6/8</b>	<b>04-06-01-21-06</b>
	<b>Tool gas tube connector 4/6</b>	<b>04-06-01-21-05</b>
	<b>Preparation needle filter</b>	<b>04-06-01-21-03</b>

## GENERAL INFORMATION

Unpack your item carefully. Inspect for damage and report such damage or missing parts to your supplier right away.

Read your Instruction Manual carefully. Take time to save time while working with your item. Make sure that every user has read and understood the Instruction Manual.

Please store the Instruction Manual in a place easily accessible to every user.

**IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !**

A so-called EURO-plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V) is standard on all of the items.

For the Continental US they feature a US-standard plug (NEMA

Pub.No.WDI.1961 ASA C 73.1 . 1961 page 8 15A 125V).

When using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one suiting your needs.

As shipped, the item features a protective ground wire. When replacing the original plug, make sure to reconnect this protective ground wire in the new plug !

## SAFETY INFORMATION

**All Safety- and Accident Prevention Regulations, as in force for laboratory work, have to be complied with.**

**When using items with built-in heaters, refer to safety information; handle with care !**

**Comply with all Safety Data Sheets published for hazardous materials. Store Safety Data Sheets in a place accessible to all operators.**

**Use eye protection and adequate garments.**

**Use extra care while working with flammable substances. Comply with Safety Data Sheets.**

**Items do not feature explosion protection.**

**Handle hot reaction block with care; danger of burns. Use protective gloves.**

**Before connecting item with your local power supply, double check with the information displayed on the data plate.**

**Turn power switch to OFF when item is not in use or before disconnecting from mains.**

**Repair work is restricted to technicians approved by Heidolph Instruments.**

**Your item requires a solid stand.**

## SET UP


**CAUTION:** keep original box for possible transit; original pack protects your item best.

## 1. Liquid-Phase Synthesis


### 1.1 Reaction Vessels

#### 1.1.1 Setup

All parts in contact with medium are made of glass or PTFE


### 1.1.2 Install head


The head is shipped as an assembly (with O-ring and cap). Head may be dismantled for cleaning. For re-assembly use tool as shown.

### 1.1.3 Built-in rotary valve

in **OPEN POSITION**, handle points towards vacuum / inert gas fitting, whereas in **CLOSED** position, handle is 90° offset.

### 1.1.4 Septa

Vessel is closed by a septum; this septum allows the injection of solvent and other substances with an injection needle during work. Septum itself is made of silicone, the inside featuring PTFE foil to ease cleaning.

### 1.1.5 Capacity of reaction vessels

For capacity of reaction vessels available for liquid-phase synthesis refer to list below:

item	P/N	total capacity	recommended capacity
reaction vessel Liquid 12	552-25030-00	50 ml	25 ml
reaction vessel Liquid 16	552-22530-00	25 ml	12 ml
reaction vessel Liquid 24	552-21030-00	10 ml	5 ml

## 1.2 Locate reaction vessels in heating- and condensation block

### 1.2.1 Remove and install reaction vessels


Handle hot reaction block with care; danger of burns. Use protective gloves.

- unthread cap from reaction vessels


- apply lateral force on reaction vessels


# E

- slip glass reaction vessel out of retainer by pushing up

- reverse above sequence of operations to install vessels
- install vessel in retainer from top end, slue heads out-of-line before

- screw head and cap on glass vessel

For perfect fixing of vessels, silicone strings are used.  
When installing silicone string, make sure to have a steady tension around the block.


right


false

When installing silicone string on reaction block Liquid 12 make sure lower string is not installed like shown on picture B.

### 1.2.2 Vacuum-/ protective gas fitting with distributor (plenum)


The vacuum- / protective gas tube connector is a snap-in type. Use tool with new tube to introduce tube in vacuum plenum.


# E

Test tube head is installed without the need for tools.

### 1.3 Connect vacuum- / inert gas lines


#### 1.3.1 Connect vacuum- / inert gas lines


The vacuum- / inert gas fitting accepts the # 8/6 PTFE tube shipped with the item. Slip union ring over PTFE tube, install clamping ring, slip hose on nipple and pull union ring tight.

#### 1.3.2 3-way valve, vacuum- / inert gas (standard equipment)

### 1.3.3 Valve settings, vacuum-/ inert gas


### 1.4 Less reaction vessels

In case of less reactions than max. capacity, leave vessels not used in reaction head and close head valve as described in Chapter 1.1.3.

## 2. Solid-phase synthesis


### 2.1 Reaction vessels

#### 2.1.1 Setup:

All parts in contact with medium are made of PFA or PTFE

**2.1.2 Install and disassemble head / base**


Slip shroud tube with clamping ring (chamfer located as shown in fig. C) over head or base (A) and seat in place (B)

Seat clamping ring against shoulder (D).

Reverse above sequence of operations for disassembly.


Head and base installed may be rotated against each other; arrange hose fittings as shown on picture in Chapter 2.1.1.

**2.1.3 Base with built-in filter**

Seat filter in base by hand. Use preparation needle coming with the item to remove filter. Mesh size is 25µm.


**2.1.4 Septa**


Reaction vessel is closed by a septum; this septum allows the injection of solvent and other substances with an injection needle during work. Septum itself is made of silicone, the inside featuring PTFE foil to ease cleaning.

### 2.1.5 Tube fitting, head and base


Slip tube fitting over PTFE tube. Slip tapered end of clamping ring over tube (refer to picture), clamp ring needs to be flush with tube end. Screw fitting with tube and clamp ring into threaded hole and pull tight by hand.

**Pull tight carefully until reaching proper vacuum seal: do not squeeze tube.**

### 2.1.6 Open and close valve of reaction vessel


## 2.1.7 Capacity of reaction vessels

Capacity of reaction vessels available for solid-phase syntheses are listed below:

item	P/N	total capacity	recommended capacity
reaction vessel Solid 16	552-14230-00	42 ml	20 ml
reaction vessel Solid 20	552-12530-00	25 ml	12 ml
reaction vessel Solid 24	552-10830-00	12 ml	5 ml

## 2.2 Locate reaction vessels in heating- and condensation block

### 2.2.1 Install vessels

Seat reaction vessels in heating -/ condensation block as shown. Exert light pressure to seat vessels. Elasticity of the shroud positively secures vessels. Vertical alignment: clamping ring is in touch with heating block.

**Handle hot reaction block with care; danger of burns. Use protective gloves.**


## 2.3 Connect vacuum / inert gas with vessels / plenum


Vacuum- / inert gas tube couplings are screw- type as shown in Chapter 2.1.5. Assignment of vessels to a particular plenum port is discretionary. However, we recommend using as short as possible tube lengths, no matter on what level you are working.

## 2.4 Connect vacuum- / inert gas to plenum ports


Vacuum- / inert gas fitting is designed for the 8/6 PTFE tube coming with the item. Slip union ring over PTFE tube, install sealing ring, plug tube over nipple, pull union ring tight.



### 2.4.1 3-way valve, vacuum-/ inert gas (standard equipment)

### 2.4.2 Valve settings, vacuum- / inert gas

## 2.5 Connect waste cube, solid-phase synthesis


Use PTFE hose (1/16") coming with the item to connect vessel base with waste cube. Make sure not to connect items with different identification markings. Install union ring and clamping ring on waste cube as shown.

Protruding length allows filtrate to be routed to its particular collecting vessel.

Waste cube vacuum supply features a tube fitting.


### 2.5.1 Convert waste cube for deviating number of vessels


To exchange waste cube adapter, loosen 3 ea. screws with hex key and remove adapter. Install new adapter, make sure seal is located between adapter and housing. Pull 3 ea. screws tight.

### 2.5.2 Schematic, arrangement of waste cube tubes


### 2.6 Less reaction vessels

Reaction block may be loaded with less reaction vessels. We however recommend not to drop below 50 % of max. capacity, in order to avoid a severe unbalance situation. Before running on vacuum or inert gas, close open ports of plenum and waste cube; to this end use plugs from the accessories list.

### 3. Coolant fitting, liquid- or solid-phase synthesis

#### 3.1 Coolant connector, condensation zone

Use high-flexible tube (6 - 8 mm ID). Secure with tube clamps.

**Do not impede shaking motion (orbit) of the reaction block by coolant tubes.**

#### 3.2 Coolant fitting, conditioning zone for reactions below room temperature.

Use high-flexible tube (6 - 8 mm ID). Insulate coolant tube when working below freezing point. Secure with tube clamps.

**Do not impede shaking motion (orbit) of the reaction block by coolant tubes**

### 4. Multiple-Evaporator

#### 4.1 Reaction vessels

refer to Chapter 1.1 liquid-phase synthesis

#### 4.2 Locate reaction vessels in heater block


refer to Chapter 1.2 liquid-phase synthesis

### 4.3 Connect vacuum

refer to Chapter 1.3 liquid-phase synthesis

### 4.4 Set-up Synthesis 1 glassware kit for Multiple-Evaporator


#### 4.4.1 Set-up stand

Screw rod in baseplate. Secure cooler holder with cross sleeve.

#### 4.4.2 Install cooler

Install cooler as shown; check for proper orientation of vacuum and coolant fittings.

#### 4.4.3 Connect vacuum and coolant

Slip special union ring over PTFE vacuum hose as shown (this special union ring is not identical to coolant caps !). Slip sealing ring in place and plug hose in vacuum fitting. Fit special union ring and pull hand tight.


#### 4.4.4 Install receiving flask

Secure receiving flask on cooler with special clamp.

#### 4.5 Schematic, Multiple-Evaporator without vacuum control

Set-up vacuum lines as shown; for more details refer to Chapter 1.3.1 and 4.4.3.

#### 4.6 Schematic, Multi-Evaporator with vacuum control

Set-up vacuum lines as shown; for more details refer to Chapter 1.3.1 and 4.4.3.

**The following components are recommended for vacuum control:**


VAC control automatic RS vacuum controller


Wouff bottle

vacuum valve

The **VAC control automatic RS** vacuum controller goes on the top end of the stand used for Synthesis 1 glassware kit.


Vacuum valve and Wouff bottle are installed on the base of the Synthesis 1 glassware kit stand (refer to Chapter 4.4).

## 5. Electric wiring

The item's rear wall bears all terminals needed.

- 1 mains (refer to voltage specification on data plate)
- 2 RS 232 serial interface
- 3 serial interface for VAC control automatic RS vacuum controller
- 4 external temperature sensor, zone 1
- 5 external temperature sensor, zone 2
- 6 external temperature sensor, zone 4
- 7 external temperature sensor, zone 3
- 8 potential-free contact to discontinue shaking action

## 6. Install guard hood

### 6.1 Guard hood


This hood is optional and protects the operator from touching the hot reaction block. Hood folds away with handle.

Handle reaction vessels with care ! Not all of the hot zones can be screened with fairings or hoods. Danger of burns !

## 6.2 Installation

Two articulated supports are installed at the rear end of the item with 2 ea. Allen bolts (1). Do not forget to install washers (2) under bolt heads. Two each adjusting screws (3) per support provide for vertical displacement of the hood as to fit different heights of reaction blocks. To this end loosen two Allen bolts (3) of both supports by abt. 1 revolution, displace as required and pull tight in position. Hex keys come with the item.


## 7. Exchange reaction blocks

Reaction blocks consist of the following components:

- 4 ea. conditioner blocks (1 ea. for every heater zone)
- condensation block
- set of reaction vessels (number depends on type of reaction block)
- vacuum- / inert gas plenum
- hoses from vacuum- /inert gas plenum to reaction vessels (number depends on type of reaction block)

First uninstall reaction block:


remove reaction vessels from reaction block


remove 2 ea. Allen bolts securing vacuum plenum on top of condensation block (hex key comes with the item), remove plenum.

# E

- remove remaining 2 screws securing condensation block; use hex key; remove condensation block

- unthread screws from 4 ea. conditioning blocks, use hex key coming with the item; remove conditioning blocks.

- remove conditioning block

- reverse above sequence of operations to install new reaction block.

## OPERATION

### 8. General


Read safety information and general notes before first operation !

Before branching item with mains supply, double check specifications on data plate. Use national standard plug. Make reference to **GENERAL NOTES** Chapter.

### 9. Controls, brief description

#### 9.1 Power switch

Item is turned **ON / OFF** with **power switch**. The item being ON, button is lighted green.

#### 9.2 Start / Stop shaking action

Shaking action is turned ON / OFF with **Start Stop Speed** key. In ON mode, a green LED turns on.

## 9.3 Start / Stop heater


Heater is turned ON / OFF with **Start Stop Heater** key. In ON mode, a green LED turns on.

## 9.4 Set / edit parameters

**Value / Speed** rotary knob sets / edits parameters. "**act** " displayed, shaking frequency is changed; **set** displayed, edits a parameter marked with an arrow key. Rotate knob clockwise to increase, counter clockwise to decrease.

**Values so set are saved automatically.**

## 9.5 Set or display

**Act / Set** key is used to select between display mode and set mode. Display mode is marked **act** in display, set mode displays **set** (set value)

## 9.6 Select value

A value to be edited in "**Set**" mode is marked by **cursor**. **Cursor** is moved with two **Select** keys.


In "**Set**" mode, **cursor** is moved with **Select** ◀ key from zone 1 to zone 2-3-4-Timer-Speed etc.;

and vice versa with **Select** ▶ key.

In "**Set Ramp**" mode, **cursor** is moved with **Select** ◀ key from No to zone – Temp Step 1 - Time Step 2 –


Temp Step 2 etc.,

and in opposite direction with **Select** ▶. key.

## 9.7 Set temperature gradient

**Set Ramp** key is used to set one or more temperature gradients


## 9.8 Activate temperature gradients

**Ramp Active** key activates a certain temperature gradient.

## 9.9 Identic

**Identic** key harmonizes 4 ea. temperature zones (all four are set to same temperature), such value being transmitted from zone 1 to the other zones. When programming temperature gradients, **Identic** key sets temperature gradient of all of the 4 heater zones.

## 10. Setting work in detail

### 10.1 Set shaking frequency

Display showing actual values ("act" displayed) and shaking is in action, rotation of **Value / Speed** knob will change shaking frequency.

Here is another option to change shaking frequency:

Select "set" with **Act / Set** key; use **Select** keys to set **cursor** to **Speed** in display and set value by rotating **Value / Speed** button. New value appears in display.

### 10.2 Set temperature (4 heater zones with different temperatures)

#### 10.2.1 Setting

Select "set" mode with **Act / Set** key; **cursor** automatically jumps to **zone 1**. Set temperature with **Value / Speed** rotary knob. New value appears in display. To set nominal values for **zone 2, 3 and 4**, set **cursor** on zone with **Select** keys and proceed as described. New values are saved right away. Return to "act" mode with **Act / Set** key.

#### 10.2.2 Temperature setting, limiting values

Temperature difference between 2 neighbored zones must not exceed 25°C. Since the four heater zones feature annular arrangement this also applies to zones at the beginning and at the end of a row.

### 10.3 Set temperature (4 heater zones with same temperature)

Select "set" with **Act / Set** key; **cursor** automatically jumps to **zone 1**. Set temperature with **Value / Speed** rotary knob. New value appears in display. To select the same nominal values for **zone 2, 3 and 4**, hit **Identic** key. Display shows same temperature for all of the zones. New values are saved automatically. To return to "act" mode, hit **Act / Set** key.

### 10.4 Set timer

A timer function allows to disconnect the heater circuit, and, after a certain and freely programmable delay, shaking action also. Timer start may be a function of activating heater circuit or activating shaker.

Select "set" mode with **Act / Set** key; move **cursor** to **Timer** with **Select** keys. Select time in minutes with **Value / Speed** rotary knob. New value appears in display.

Heater circuit is disconnected as programmed, shaking action after a certain delay, freely programmable in the **Setup** mode (refer to Chapter 12).

You may also activate **Timer** while synthesis is in process; in this case, time remaining is set.

**Timer** being active is confirmed by a moving arrow beside time displayed.

**Timer** not being activated, display (in display mode, Timer mode) shows time elapsed since start of heater mode or start of shaking action. **Timer** start may be triggered by heater start or shaking action start. Select in **Setup** mode (refer to Chapter 12). Return to "act" mode with **Act / Set** key.

## 10.5 Set temperature gradient

### 10.5.1 Set values

A temperature gradient for an individual heater zone may be programmed in up to 9 steps. Temperature changes from initial step to next step within a period of time selected before. Times to be set accumulated progress steps and not sectional values. (refer to example)

A total of 4 different temperature gradients may be saved.

Hit **Set Ramp** key. A temperature gradient setting menu appears in display. **Cursor** automatically jumps to No X. This is the file name of the previous parameter set. To select another of the 4 memories, or overwrite a memory, select this memory with the **Value** rotating knob.

Move **cursor** to zone wanted with **Select** keys. Select zone to be edited with **Value** rotary knob. To program all of the 4 zones with the same temperature gradient, hit **Identic** key. Display changes from zone X to "Ident".

At this point set cursor on first temperature in table and set value by rotating Value knob. ( Start-time is "0" all the time).

**Select** ◀ **key** moves cursor from field to field. **Select** ▶ **key** takes you to next zone automatically when programming is complete.

Set Ramp	No X	Zone X							
Step	1	2	3	4	5	6	7	8	9
Time	0	---	---	---	---	---	---	---	---
Temp	▶0	0	0	0	0	0	0	0	0

Click on more values with cursor and set value with rotary knob. Refer to physical limiting values for heating and cooling.

**Time value of a step not being „---“, , number zero will be counted as 0°C !**

Put a straight line through two values of a sequence (refer to graph)

max. heating rate: 3°C per minute

max. cooling rate in air: 2°C per minute

Examples of a ramp so programmed:


Set Ramp	No X	Zone X							
Step	1	2	3	4	5	6	7	8	9
Time	0	10	100	120	200	250	---	---	---
Temp	20	40	40	65	65	▶25	---	---	---

Temperature value of the last step programmed will be maintained until disconnecting the item from power supply. Item accepts up to 9 steps, less of course also. To save room in display, only 3 steps are displayed the same time; to have a look at more steps, use cursor scroll function (Select key). e.g. after steps 1-2-3, you may scroll to steps 2-3-4 or 4-5-6, etc. To return, hit opposite-sense **Select** key. Not needing to program all of the 9 steps, hit **Set Ramp** key to directly reach the next zone's programming function. No temperature changes wanted for certain zones, just enter nominal temperature for the zone concerned in step 1.

Hit **Set Ramp** key once again to exit gradient setting.

Having set a non-plausible time value, you may not exit temperature gradient setting mode with **Set Ramp** key. Display shows **“time invalid”** and cursor marks non-plausible value that needs to be changed. Time values from step 1 to step 9 need to be ascending.

### 10.5.2 Delete values

For deleting a cursor-marked time value, rotate "Value" button CCW until „---“ is displayed. Adhering temperature value automatically is set to „0“.

### 10.5.3 Limiting temperatures

Within a temperature range of 20°C and 160°C, max. temperature rise / drop, in absence of optional forced cooling, makes reference to the table below (for physical reasons):

temperature range	max. temperature rise	max. temperature drop
20 – 80 °C	3°C / min	1°C / min
81 – 130 °C	2 °C / min	2°C / min
131 – max °C	1°C / min	3°C / min

Deviating temperature rise values cannot be achieved within the time interval specified; the same of course is true for temperature drop.

## 10.6 Activate temperature gradient

Activate temperature gradient by hitting the **Ramp active** key. This mode is confirmed by the **Ramp active** LED. A temperature gradient set as described in

Chapter 10.5.1 or recalled from memory will be run when heater is started. Actual temperature of individual zones may be monitored in the “**Act**” mode. In “**Set**” mode, the actual and calculated “**Set**” value can be monitored.

When ramp is activated, settings can be viewed by pressing **Set Ramp**. values can not be changed in this mode!


## 10.7 Turn shaking action ON

**Start Stop Speed** key turns shaking action ON/OFF. Green LED is lighting shaking action being ON.

## 10.8 Turn heater ON

**Start Stop Heater** turns the heater ON/OFF. Green LED is lighting heater being ON.

Current dissipation in one of the heating zones is displayed by a symbol beside the temperature in “**Act**” mode.

## 11. External temperature sensor

An external temperature sensor measures or controls temperature in a reaction vessel.

One external temperature may be connected with one heating zone. This sensor is introduced through septa. To this end first pre-punch septa with a thick injection needle (refer to Chapter 1.1.4 and 2.1.4 Septa) Introduce sensor through pre-punched hole careful.

**IMPORTANT: double-check proper assignment of temperature sensors to heating zones.**

Re-adjust vertical position of hood as required (refer to Chapter 6, Set-up)

It is important that even while shaking, 20 mm of the front section of the sensor is wetted by the reaction fluid.



## 11.1 Temperatures below solvent's boil point

While working at temperatures below solvent's boil point, temperature sensor may be used to regulate temperature in the reaction vessel.

**Important:** to this end set "control" in Setup mode (refer to Chapter 12)

An external temperature sensor being connected, display in the "Act" mode marks every individual zone with this symbol. This symbol appearing, temperature of this zone is regulated by the external temperature sensor. In this case temperature of the heating block will be higher than actual temperature measured, in order to compensate for physical temperature decrease.


## 11.2 Temperatures above solvent's boil point.

While working at temperatures above solvent's boil point, temperature sensor may be used to measure temperature in the reaction vessel only.

**IMPORTANT:** to this end select "report" in the setup mode (refer to Chapter 12)

External temperature sensor being connected, display in the "Act" mode is marked with this symbol for every single zone. This symbol appearing means that temperature in the reaction vessel of this zone is displayed.

In "Act" mode, temperature of the conditioner block may be displayed either. To this end hit one of the **Select** keys. Temperature of the conditioner block is marked with this symbol.

## 11.3 In the use as multi evaporator

When evaporating, use temperature sensors in report-mode only (see chap. 11.2). Else sample could overheat.

## 12. Setup (timer and temperature sensor)

In the setup mode the following parameters may be set:

- Start Timer by turning heater or shaking action ON.
- Such time being elapsed, heating action and shaking action are stopped the same time; as another option, heater stops first, followed by shaking motion after a certain delay.
- Sensor

### How to proceed:

Hold **Set/Act** key depressed and turn power switch to **ON**.  
The following information appears in display:

```
SETUP
Start timer
▶ Heater
Speed
```

Move cursor to place required with **Select** keys, thereafter hit **Set/Act** key  
The following information appears in display:  
At this point you may delay time for shaking action (after heater has been

```
SETUP
Stop timer
Speed:
▶ 000 min time delay
```

disconnected). Such time being set by rotating the **Value** knob between 0 and 999 minutes.

thereafter hit **Set** key.  
The following information appears in display:


```

SETUP
Sensor
▶control
report
  
```

Put **cursor** to position wanted with **Select** keys, hit **Set/Act** key. (refer to Chapter 11.1 and 11.2)

At this point set-up is complete, display returns to standard screen.

### 13. Filter action, solid-phase synthesis

- Make sure vessels are aerated via the vacuum / inert gas port.

- Open waste cube and load with beaker. Mechanically close waste cube.
- Open base valves in reaction vessels (refer to Chapter 2.1.6). Evacuate waste cube (e.g. start vacuum pump).
- Contents of reaction vessel now is drawn-off via filters and collected in beaker.
- Connect waste cube with atmospheric pressure (turn vacuum pump OFF)
- Close valves in reaction vessels (refer to Chapter 2.1.6).
- After having separated product from resin, remove beaker from waste cube and fit test tube rack in waste cube. To this end squeeze grip zones of rack and load rack, still holding the marked areas squeezed. Mechanically close waste cube. From this point on proceed as described.

## 14. Use condensation zone

Running above boiling point of the reaction fluid, vapor may be separated in a condensation zone. To this end connect condensation zone to coolant circuit or a recirculating cooler as described in Chapter 3.1.

## 15. Use heating zones with recirculating cooler (low temperature)

Working with reactions at temperatures below ambient, use an external recirculating cooler for the conditioning zone. To this end connect recirculating cooler to conditioning zone as described in Chapter 3.2 before starting work.

When working at temperatures below freezing, expect significant accumulation of ice on the reaction block. Ice impedes view and makes item slippery (handle with care) !

**Never attempt to work with item at temperatures below freezing in absence of protective gas atmosphere. Locate item in drip pan and fill with protective gas.**


Temperatures below freezing also may be run by the different heating zones. Coolest zone is close to the coolant entry port (e.g. zone 4); rise temperatures in the following sequence: 4 - 1 - 2 - 3. (Using coolant entry port in zone 3, use 3 - 2 - 1 - 4 sequence). Different temperatures are reached by heating "against" coolant action. To this end select temperatures and start heater as described in Chapter 10.2 and 10.3.

**Having worked at low temperatures and switching over to high temperatures requires to drain all of the coolant from its passage.**

## 16. Access to rear reaction vessels

To increase view on rear reaction vessels in a block, reaction block slues through 30 degrees (both sides). To this end push or pull lever on the right of the reaction block). Return reaction block to central position before starting shaking. In central position, lever drops with an audible click.

## 17. Interface RS 232

### Interface commands

#### a. Query

SZn? with n=1...4                      nominal temperature, zone n  
*answer: SZn=xxx*

SD?    nominal rotational speed  
*answer: SD=xxxx*

ST?    nominal time, timer elapse  
*answer: ST=xxx*

BT?    running down, timer  
*answer: BT=xxx*

AZn? with n=1...4                      actual temperature, zone n  
*answer: AZn=xxx*

AEn? with n=1...4                      actual value, external sensor, zone n  
*answer: AEn=xxx*

AD?    actual rotational speed  
*answer: AD=xxxx*

SB?    operational status  
*answer: SB=x*

Answer	heater	speed	Ramp Active	Identic
Bit	b0	b1	b2	b3
SB=0	off	off	off	off
SB=1	on	off	off	off
SB=2	off	on	off	off
SB=3	on	on	off	off
SB=4	off	off	on	off
SB=5	on	off	on	off
SB=6	off	on	on	off
SB=7	on	on	on	off
SB=8	off	off	off	on
SB=9	on	off	off	on
SB=10	off	on	off	on
SB=11	on	on	off	on
SB=12	off	off	on	on
SB=13	on	off	on	on
SB=14	off	on	on	on
SB=15	on	on	on	on

FN?    error number  
*answer: FN=x*

RZ? ramp values, time  
*answer:* RZ(x,y)=aaa,bbb,ccc,ddd,eee,fff,ggg,hhh,iii  
 Note: ramp number / zone need to be set before  
 x=ramp number,y=ramp zone

RT? ramp values, temperature  
*answer:* RT(x,y)=aaa,bbb,ccc,ddd,eee,fff,ggg,hhh,iii  
 Note: ramp number / zone need to be set before  
 x=ramp number,y=ramp zone

AP? actual pressure  
*answer:* AP=xxxx

**b. Entry**

SZn=xxxx! with n=1...4 nominal temperature, zone n  
*answer:* SZn=xxx

SD=xxxx! nominal rotational speed  
*answer:* SD=xxxx

ST=xxx! nominal time, timer elapse  
*answer:* ST=xxx

SB=xxx! operational status  
*answer:* SB=xxx

Entry	heater	speed	ramp active	Identic
Bit	b0	b1	b2	b3
SB=000!	off	off	off	off
SB=001!	on	off	off	off
SB=002!	off	on	off	off
SB=003!	on	on	off	off
SB=004!	off	off	on	off
SB=005!	on	off	on	off
SB=006!	off	on	on	off
SB=007!	on	on	on	off
SB=008!	off	off	off	on
SB=009!	on	off	off	on
SB=010!	off	on	off	on
SB=011!	on	on	off	on
SB=012!	off	off	on	on
SB=013!	on	off	on	on
SB=014!	off	on	on	on
SB=015!	on	on	on	on

RN=x! ramp number  
*answer:* RN=x

RR=x! ramp zone  
*answer:* RR=x  
 x=1...4 zone 1...4  
 x=5 Identic

# E

RZ= aaa,bbb,ccc,ddd,eee,fff,ggg,hhh,iii  
ramp values, time

*answer:* RZ(x,y)=aaa,bbb,ccc,ddd,eee,fff,ggg,hhh,iii

Note: ramp number / zone need to be set before  
x=ramp number,y=ramp zone

RT= aaa,bbb,ccc,ddd,eee,fff,ggg,hhh,iii  
ramp values, temperature

*answer:* RT(x,y)=aaa,bbb,ccc,ddd,eee,fff,ggg,hhh,iii

Note: ramp number / zone need to be set before  
x=ramp number,y=ramp zone

**Interface parameters:** 9600 Baud, no parity, 8 bit, 1 stop bit.

- RS 232 interface being active, all controls and displays of the item remain active and serviceable.

## CLEANING & SERVICING

**Cleaning:** wipe housing surfaces clean with a damp cloth (add some sort of mild liquid soap).

### Note

To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based detergents, abrasive substances, ammonia, rags or cleaning agents containing metal particles.

The item is maintenance-free. Repair work is limited to technicians authorized or appointed by Heidolph Instruments. Please call your local Heidolph Instruments Dealer or a Heidolph Instruments Field Representative

## UNINSTALL, FORWARD & STORE

### Disassembly

1. Turn item OFF and disconnect mains plug.
2. Remove all of the hardware arranged around the shaker to ease disassembly.
3. Unload all reaction vessel from reaction block, uninstall optional equipment.

### Forward & Store

1. We recommend to store the item and its components in its original box, or a similar container that offers adequate protection against damage in transit. Tape the box securely.
2. Store the item in a dry place.

**CAUTION: Do not jolt or shake the item during transport.**

**CAUTION: Use original box to avoid damage in transit.**

## DISPOSAL

For disposal, please comply with your local or national regulations.  
Split by metal, plastic, etc.  
Packing material to be treated as described above (material split).

## TROUBLESHOOTING

Work on electric, electronic and cryogenic components is limited to qualified personnel.

### **Error messages in display:**

One of the following error messages appearing in display, the item is disconnected automatically for reasons of safety.

**ERROR**  
Sensor intern zone X

Temperature sensor in heating zone X (1 / 2 / 3 or 4) broken.  
Replace temperature sensor (Heidolph After Sales Service).

**ERROR**  
Overheat protection

Overtemperature circuit breaker dropped. Ship item to Heidolph After Sales Service for repair.

**ERROR**  
System

Defect in electronic circuit. Ship item to Heidolph After Sales Service for repair.

# ERROR

## Sensor extern zone X

Extern temperature sensor in zone X (heating zone 1 / 2 / 3 or 4) broken.  
Replace temperature sensor.

### SPECIFICATIONS

Model	<i>Synthesis 1</i> Solid Phase Synthesis	<i>Synthesis 1</i> Liquid Phase Synthesis	<i>Synthesis 1</i> Multiple Evaporator
P/N	Synthesis 1 Solid 16: 552-14200-00 Synthesis 1 Solid 20: 552-12500-00 Synthesis 1 Solid 24: 552-10800-00	Synthesis 1 Liquid 12: 552-25000-00 Synthesis 1 Liquid 16: 552-22500-00 Synthesis 1 Liquid 24: 552-21000-00	Multiple Evaporator 12: 552-35000-00 Multiple Evaporator 16: 552-32500-00 Multiple Evaporator 24: 552-31000-00
Footprint (wxd) [mm]	310 x 500	310 x 500	310 x 500
Material of reaction vessels	transparent PFA & PTFE	glass, PTFE caps	glass, PTFE caps
Vibration frequency [rpm]	100 – 1000	100 – 1000	100 – 1000
Orbit [mm]	3	3	3
Drive	brushless, electronic D. C. motor	brushless, electronic D. C. motor	brushless, electronic D. C. motor
Heating	4 different heating zones with individual heat control and temperature display	4 different heating zones with individual heat control and temperature display	4 different heating zones with individual heat control and temperature display
Heating power [W]	approx. 1000	approx. 1000	approx. 1000
Temperature range / heating [°C]	up to 160	up to 160	up to 160
Temperature range / cooling	down to - 80 by external chiller	down to - 80 by external chiller	down to - 80 by external chiller
Temperature accuracy [°C]	1	1	1
Temperature display	digital	digital	digital
Temperature ramping	4 heating zones with individual heat control allow for programming of automated temp. ramping with up to 9 steps	4 heating zones with individual heat control allow for programming of automated temp. ramping with up to 9 steps	4 heating zones with individual heat control allow for programming of automated temp. ramping with up to 9 steps
Number of reaction vessels [max.]	16 reaction vessels of 42 ml each 20 reaction vessels of 25 ml each 24 reaction vessels of 8 ml each	12 reaction vessels of 50 ml each 16 reaction vessels of 25 ml each 24 reaction vessels of 10 ml each	12 reaction vessels of 50 ml each 16 reaction vessels of 25 ml each 24 reaction vessels of 10 ml each
Connection for inert gas atmosphere / vacuum	via distributor block	via distributor block	via distributor block
Adding of solvent	via injection needles into septa	via injection needles into septa	via injection needles into septa
Capacity of Waste Cube	approx. 1000 ml	no Waste Cube	no Waste Cube
Interface RS 232	for monitoring and controlling temperature, vibration frequency and pressure*	for monitoring and controlling temperature, vibration frequency and pressure*	for monitoring and controlling temperature, vibration frequency and pressure*
Supply voltage	230 V 50/60 Hz 115 V 50/60 Hz 100 V 50/60 Hz	230 V 50/60 Hz 115 V 50/60 Hz 100 V 50/60 Hz	230 V 50/60 Hz 115 V 50/60 Hz 100 V 50/60 Hz
Operating conditions [°C]	0 - 40 at 80 % relative humidity	0 - 40 at 80 % relative humidity	0 - 40 at 80 % relative humidity
Protection class	IP 20	IP 20	IP 20
Weight [kg]	27,5	27,5	28,5
Timer [min]	0 - 999	0 - 999	0 - 999
Over temperature protection	electronic	electronic	electronic
Supply power [W]	1120	1120	1120

\* pressure via optional vacuum controller

## WARRANTY, LIABILITY & COPYRIGHT

### Warranty

Heidolph Instruments warrants that the present product shall be free from defects in material (except wear parts) and workmanship for 3 years from the date shipped off the manufacturer's warehouse.

Transit damage is excluded from this warranty.

To obtain such warranty service, contact Heidolph Instruments (phone: +49 – 9122 - 9920-67) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge. Misuse, abuse, neglect or improper installation are not covered by this warranty promise.

Alterations to the present warranty promise need Heidolph Instruments' consent in writing.

### Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

### Copyright

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

## QUESTIONS / REPAIR WORK

If any aspect of installation, operation or maintenance remains unanswered in the present Manual, please contact us at the following address:

For equipment repairs, please call Heidolph Instruments (phone: +49 – 9122 - 9920-67) or your local authorized Heidolph Instruments Dealer.

#### Note:

You will receive approval for sending your defective item to the following address:

#### Note:

If you are based in the United States of America, please contact Heidolph US:

**Heidolph Instruments GmbH & Co. KG**  
**Lab Equipment Sales**  
**Walpersdorfer Str. 12**  
**D-91126 Schwabach / Germany**  
**Tel.: +49 – 9122 - 9920-68**  
**Fax: +49 – 9122 - 9920-65**  
**E-Mail: [sales@heidolph.de](mailto:sales@heidolph.de)**

**Heidolph Instruments, LLC**  
**Lab Equipment Sales**  
**2615 River Rd.**  
**Cinnaminson, NJ 08077**  
**Phone: 856-829-6160**  
**Fax: 856-829-7639**  
**E-Mail: [heidolph@snip.net](mailto:heidolph@snip.net)**

**Safety Information**

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take adequate *protective* action not to harm our parts receiving and service personnel
- *mark* the pack IAW Hazardous Materials Act

**CE-DECLARATION OF CONFORMITY**

We herewith declare that the present product complies with the following standards and harmonized documents:

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Respuesta	Calefacción	Speed	RampActive	Identic
Bit	b0	b1	b2	b3
SB=0	off	off	off	off
SB=1	on	off	off	off
SB=2	off	on	off	off
SB=3	on	on	off	off
SB=4	off	off	on	off
SB=5	on	off	on	off
SB=6	off	on	on	off
SB=7	on	on	on	off
SB=8	off	off	off	on
SB=9	on	off	off	on
SB=10	off	on	off	on
SB=11	on	on	off	on
SB=12	off	off	on	on
SB=13	on	off	on	on
SB=14	off	on	on	on
SB=15	on	on	on	on

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Comando	Calefacción	Speed	RampActive	Identic
Bit	b0	b1	b2	b3
SB=000!	off	off	off	off
SB=001!	on	off	off	off
SB=002!	off	on	off	off
SB=003!	on	on	off	off
SB=004!	off	off	on	off
SB=005!	on	off	on	off
SB=006!	off	on	on	off
SB=007!	on	on	on	off
SB=008!	off	off	off	on
SB=009!	on	off	off	on
SB=010!	off	on	off	on
SB=011!	on	on	off	on
SB=012!	off	off	on	on
SB=013!	on	off	on	on
SB=014!	off	on	on	on
SB=015!	on	on	on	on

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